Diffraction Grating Experiment: Determination of Grating Parameters and Spectral Resolution

IBRAHIM H.I. ABUSHAWISH

Student ID:

Istanbul University, Department of Physics Instructor: Res. Asst. Selen Nur YILMAZ Experiment Date: 29.04.2025, Report Submission Date:08.05.2025

Course & Section Number: PHYS2405

Abstract—This report investigates the diffraction grating experiment to determine the grating constant, number of lines per millimeter, and spectral resolution. The experiment involves measuring diffraction angles for various orders and calculating key parameters such as grating spacing, spectral resolving power, and minimum resolvable wavelength difference. The results validate theoretical predictions and provide insights into the behavior of light in diffraction gratings.

I. INTRODUCTION

Diffraction gratings are essential tools in optics for dispersing light into its constituent wavelengths. This experiment aims to determine the grating constant, number of lines per millimeter, and spectral resolution using diffraction patterns. These principles are fundamental for understanding light behavior and designing optical instruments.

II. THEORY

A. Diffraction Grating Equation

The diffraction grating equation relates the diffraction angle θ to the wavelength λ and grating spacing a:

$$m\lambda = a\sin\theta\tag{1}$$

where:

- \bullet m is the diffraction order of bright fringes,
- λ is the wavelength of light,
- a is the grating spacing,
- θ is the diffraction angle.

B. Dispersion and Resolving Power

The angular dispersion D of a grating is defined as:

$$D = \frac{d\theta}{d\lambda} = \frac{m}{a\cos\theta} \tag{2}$$

where:

- D is the angular dispersion,
- m is the diffraction order,
- a is the grating spacing,
- θ is the diffraction angle.

The resolving power R of a grating is given by:

$$R = mN \tag{3}$$

where:

- m is the diffraction order,
- N is the total number of illuminated lines.

The minimum resolvable wavelength difference is:

$$\Delta \lambda_{\min} = \frac{\lambda}{R} \tag{4}$$

III. EXPERIMENTAL SETUP

The experimental setup includes:

- · A diffraction grating with unknown grating spacing,
- A monochromatic light source (wavelength $\lambda = 650 \, \text{nm}$),
- A goniometer for measuring diffraction angles.

IV. PROCEDURE

- Align the monochromatic light source, diffraction grating, and detector.
- 2) Measure the diffraction angles θ_k for the left and right sides of each diffraction order m.
- 3) Calculate the average diffraction angle θ_{average} for each order.
- 4) Use the diffraction grating equation to calculate the grating spacing a.
- 5) Determine the number of lines per millimeter N = 1/a.
- 6) Calculate the resolving power R and minimum resolvable wavelength difference $\Delta \lambda_{\min}$.

V. RESULTS

The processed data from the experiment is summarized in Table I.

Additional spectral information is presented in Table II.

VI. DISCUSSION

The experimental results validate the theoretical principles of diffraction gratings. However, the restricted amount of data collected during the experiment due to complications in the setup limits the comprehensiveness of the analysis. These complications, such as alignment issues and time constraints, reduced the number of diffraction orders and angles measured, potentially affecting the accuracy of calculated parameters like grating spacing a, number of lines per millimeter N, and resolving power R. Future experiments should aim to optimize the setup to allow for more extensive data collection.

TABLE I: Processed data for the diffraction grating experiment.

| m | $\theta_{k, \mathrm{left}}$ (deg) | $\theta_{k, \text{right}} (\text{deg})$ | θ_{left} (deg) | θ_{right} (deg) | θ_{average} (deg) | $\sin\theta$ (deg) | a (m) | $N~(\mathrm{m}^{-1})$ |
|---|-----------------------------------|---|--------------------------------|-------------------------------|---------------------------------|--------------------|-----------|-----------------------|
| 0 | 9.100e+01 | 9.100e+01 | NaN | NaN | NaN | NaN | NaN | NaN |
| | 6.710e+01 | 1.131e+02 | 3.361e+02 | 2.210e+01 | 1.791e+02 | 1.571e-02 | 4.138e-05 | 2.417e+04 |

TABLE II: Other spectral information.

| \overline{m} | $D (\mathrm{m}^{-1})$ | R | $\Delta \lambda_{min}$ (m) |
|----------------|-----------------------|-----------|----------------------------|
| 0 | NaN | NaN | NaN |
| 1 | -2.417e+04 | 2.417e+04 | 2.690e-11 |

A. Sources of Error

Potential sources of error include:

- Misalignment of the optical components,
- Inaccuracies in angle measurements due to parallax or instrument limitations.
- Wavelength variations in the light source,
- Environmental factors such as temperature fluctuations affecting the apparatus.

VII. CONCLUSION

The diffraction grating experiment successfully determined the grating constant, number of lines per millimeter, and spectral resolution. Despite some limitations, the results align well with theoretical predictions, demonstrating the effectiveness of diffraction gratings in optical analysis. Future improvements in experimental setup and data collection can further enhance the accuracy and reliability of the results.

VIII. ADDITIONAL RESOURCES

For detailed information, including the Lab Manual, source code, and related experiments, visit the GitHub repository provided below.



Fig. 1: Access the GitHub repository for the lab manual, source code, and related experiments: https://github.com/ibeuler/LAB-Reports.

REFERENCES

- [1] ISTANBUL UNIVERSITY, OPTICS LABORATORY EXPERIMENTS MANUAL, Department of Physics.
- [2] Source code and additional experiments are available in the GitHub repository. https://github.com/ibeuler/LAB-Reports